

CLEMENT CLARKE INTERNATIONAL PRESS RELEASE

Clement Clarke International Ltd. • Edinburgh Way • Harlow • Essex CM20 2TT • Phone 01279 414969 • Fax 01279 456305

FOR IMMEDIATE RELEASE

Date: 15th September 2005

Phone: 01279 414969

Fax: 01279 456305

CLEMENT CLARKE INTERNATIONAL LAUNCHES SPACE CHAMBER

Harlow, Essex -- Clement Clarke International Ltd, respiratory manufacturer of well-established, trusted brands such as Mini-Wright, Medix and Lifecare, will launch the new Space Chamber, spacer device, at the European Respiratory Society Congress in Copenhagen, Denmark, September 17th – 21st 2005.

The Space Chamber is a spacer device which is designed to make it easier and more efficient to use an inhaler. By providing a space between the inhaler and the mouth, the Space Chamber reduces the velocity of the aerosol and the subsequent impact on the throat. Unwanted large particles in the aerosol mist are trapped inside the Space Chamber preventing irritation of the mouth and throat. In addition it allows more time for evaporation of the propellant so that a larger proportion of the particles can be inhaled and deposited on the lungs. Using the Space Chamber, in both adults and children, more of the right size particles of medication enter deep into the lungs than when using an inhaler on its own. The Space Chamber is also available with an oxygen port for convenient delivery of aerosol medication without interrupting oxygen therapy.

The Space Chamber joins Clement Clarke's Able Spacer product in their spacer device offering. The Able Spacer is a discreet, compact small-volume spacer, with a built-in coaching signal that helps patients maintain proper inspiratory rate for maximum drug delivery. Cleverly designed to house the complete inhaler, when not in use, the universal adaptor port and protective flap accommodates all commonly used metered dose inhalers.

Clement Clarke International Ltd is dedicated to combining years of respiratory experience and the proven performance of trusted brands, with the distinctive talent for offering unparalleled support and tailoring customised solutions to enhance patient care, while improving efficiency and reducing cost. The company continues to develop proven brand product ranges to strengthen our respiratory portfolio in asthma management, COPD, oxygen therapy, cystic fibrosis and allergic rhinitis.

For more information on both the Space Chamber and Able Spacer, contact customer services on (01279) 414969.

-ENDS-

For further press information, in the first instance please contact:

Helen Powell, Clement Clarke International Ltd. Tel: 01279 456247