

# CLEMENT CLARKE INTERNATIONAL PRESS RELEASE

---

Clement Clarke International Ltd. • Edinburgh Way • Harlow • Essex CM20 2TT • Phone 01279 414969 • Fax 01279 456305

**FOR IMMEDIATE RELEASE**

**Date: 15<sup>th</sup> September 2005**

**Phone: 01279 414969**

**Fax: 01279 456305**

## **CLEMENT CLARKE INTERNATIONAL OFFERS TWIN SCALE MINI-WRIGHT FOR EASY RESULTS CONVERSION**

Harlow, Essex -- Clement Clarke International Ltd, respiratory manufacturer of well-established, trusted brands such as Mini-Wright, Medix and Lifecare, now offers a Twin Scale Mini-Wright peak flow meter for easy conversion of results.

With the introduction of the new EU 13826 peak expiratory flow standard, Clement Clarke has developed a new Twin Scale version of the standard Mini-Wright peak flow meter. The twin scale peak flow meter features both the new EU scale and the Wright scale side by side for fast, simple conversion of PEF results.

Measuring PEF between 60-800 l/min, the Twin Scale Mini-Wright displays the new EU scale, developed to improve the management of respiratory disease and adopted across Europe in 2004, in yellow with blue text. The Wright-McKerrow scale is presented in black with white numbers.

Clement Clarke International Ltd is dedicated to combining years of respiratory experience and the proven performance of trusted brands, with the distinctive talent for offering unparalleled support and tailoring customised solutions to enhance patient care, while improving efficiency and reducing cost. The company continues to develop proven brand product ranges to strengthen our respiratory portfolio in asthma management, COPD, oxygen therapy, cystic fibrosis and allergic rhinitis.

For more information on the Twin Scale Mini-Wright peak flow meter, contact customer services on (01279) 414969.

**-ENDS-**

**For further press information, in the first instance please contact:**

Helen Powell, Clement Clarke International Ltd. Tel: 01279 456247